A Novel Articulating Endoscopic Tandem and Inflatable Ovoids for the Intracavitary Brachytherapy of Cervical and Uterine Cancers

- SPECIFICATIONS
- The endoscopic dual-lumen articulating tandem is introduced into the cervical os and advanced into the uterine canal. The tandem is precisely adjusted to align with the uterine axis, ensuring optimal and reproducible inter-fraction implant geometry.
- **Insertion of Inflatable Ovoids:** A deflated dual-ovoid balloon assembly is positioned within the vaginal fornices. Once in place, the ovoids are expanded in the horizontal plane to achieve the desired spatial distribution within the pelvic cavity.
- **Insertion of Bladder and Rectal Retractors:** The bladder and rectal retractors are sequentially inserted. Their lengths are meticulously adjusted to accommodate patient-specific anatomy and subsequently secured in position.
- **Inflation of Balloons**: The ovoid and retractor balloons are inflated with sterile saline or water under controlled pressure via a one-way infusion valve. The inflation process is conducted externally, without direct internal patient manipulation.

- **Pre-sterilized, Single-Use Design:** Provided in a sterile, ready-to-use configuration to ensure aseptic conditions and eliminate the need for reprocessing.
- **High-Resolution CMOS Camera Integration:** Enables real-time visualization for precise tandem placement, minimizing the risk of uterine perforation.
- **Anatomically Adaptive Inflatable Ovoids:** Designed to conform to patient-specific anatomy, reducing trauma associated with rigid plastic ovoids and enhancing radiation dosimetry through differential loading.
- **Optional Inflatable Retractors:** Replace traditional gauze packing with a standardized approach to protect the bladder and rectum, ensuring reproducibility in applicator placement.

**Multi-Modality Imaging Compatibility:** Fully compatible with MRI, CT, and X-ray imaging for accurate treatment planning and verification.

- **Extracorporeal Inter-Fraction Geometric Optimization:** Implant adjustments are performed externally with precision comparable to laparoscopic techniques, ensuring reproducibility across treatment fractions.
- **Articulating Tandem Adaptability:** The tandem conforms to the natural anatomical angulation of both anteverted and retroverted uteri, facilitating precise tumor targeting while minimizing radiation exposure to adjacent critical structures.
- **Adjustable Implant Geometry for Conformal Therapy:** This approach allows for dynamic modification of applicator positioning, optimizing radiation dose distribution to achieve highly conformal brachytherapy.
- **Reduced Procedural Complexity and Anesthesia Requirements:** By streamlining the implantation process, this method has the potential to decrease procedural duration, thereby reducing the need for prolong sedation or general anesthesia and enhancing patient comfort.

**Cost-Effective and Consistent Brachytherapy Implementation:** The technique may eliminate the necessity for Smit-Sleeve placement and improve the reproducibility of applicator positioning, leading to greater efficiency and cost savings in clinical practice.

**Briefing Document: Novel Articulating Endoscopic Tandem and Inflatable Ovoids for Intracavitary Brachytherapy**

**Executive Summary:**
This document reviews a novel brachytherapy device and technique for treating cervical and uterine cancers. The core innovation revolves around an articulating endoscopic tandem coupled with inflatable ovoids and retractors. This system aims to improve precision, reduce trauma, enhance treatment conformity, simplify procedures, and potentially reduce costs compared to traditional methods. Key features include real-time visualization via a CMOS camera, anatomically adaptive inflatable components, compatibility with multiple imaging modalities, and the ability to adjust implant geometry externally between treatment fractions.
**Key Themes and Ideas:**

1. **Improved Precision and Reproducibility:** The system emphasizes precise placement and reproducible inter-fraction geometry. The articulating tandem is designed to align with the uterine axis, and the inflatable ovoids can be positioned accurately within the vaginal fornices.

- Quote: "*The endoscopic dual-lumen articulating tandem is introduced into the cervical os and advanced into the uterine canal. The tandem is precisely adjusted to align with the uterine axis, ensuring optimal and reproducible inter-fraction implant geometry.*"
- Quote: "*Extracorporeal Inter-Fraction Geometric Optimization: Implant adjustments are performed externally with precision comparable to laparoscopic techniques, ensuring reproducibility across treatment fractions.*"

1. **Reduced Trauma and Enhanced Patient Comfort:** Inflatable ovoids and retractors are designed to conform to patient anatomy, replacing traditional rigid components and gauze packing, potentially minimizing trauma. The streamlined implantation process also aims to decrease procedural duration, thereby reducing the need for prolonged sedation or general anesthesia.

- Quote: "*Anatomically Adaptive Inflatable Ovoids: Designed to conform to patient-specific anatomy, reducing trauma associated with rigid plastic ovoids and enhancing radiation dosimetry through differential Radiation loading.*"
- Quote: "*Optional Inflatable Retractors: Replace traditional gauze packing with a standardized approach to protect the bladder and rectum, ensuring reproducibility in applicator placement.*"
- Quote: "*Reduced Procedural Complexity and Anesthesia Requirements: By streamlining the implantation process, this method has the potential to decrease procedural duration, thereby reducing the need for prolong sedation or general anesthesia and enhancing patient comfort.*"

1. **Conformal Therapy and Optimized Dose Distribution:** The adjustable implant geometry allows for dynamic modification of applicator positioning, optimizing radiation dose distribution to achieve highly conformal brachytherapy. The articulating tandem can adapt to different uterine orientations, minimizing radiation exposure to adjacent critical structures.

- Quote: "*Adjustable Implant Geometry for Conformal Therapy: This approach allows for dynamic modification of applicator positioning, optimizing radiation dose distribution to achieve highly conformal brachytherapy.*"
- Quote: "*Articulating Tandem Adaptability: The tandem conforms to the natural anatomical angulation of both anteverted and retroverted uteri, facilitating precise tumor targeting while minimizing radiation exposure to adjacent critical structures.*"

1. **Simplified Workflow and Potential Cost Savings:** The pre-sterilized, single-use design eliminates the need for reprocessing, ensuring aseptic conditions. The technique potentially eliminates the necessity for Smit-Sleeve placement and improves the reproducibility of applicator positioning, leading to greater efficiency and cost savings.

- Quote: "*Pre-sterilized, Single-Use Design: Provided in a sterile, ready-to-use configuration to ensure aseptic conditions and eliminate the need for reprocessing.*"
- Quote: "*Cost-ENective and Consistent Brachytherapy Implementation: The technique may eliminate the necessity for Smit-Sleeve placement and improve the reproducibility of applicator positioning, leading to greater eMiciency and cost savings in clinical practice.*"

1. **Real-time Visualization and Multi-Modality Imaging Compatibility:** The integrated high-resolution CMOS camera enables real-time visualization for precise tandem placement, minimizing the risk of uterine perforation. The system is also fully compatible with MRI, CT, and X-ray imaging for accurate treatment planning and verification.

- Quote: "*High-Resolution CMOS Camera Integration: Enables real-time visualization for precise tandem placement, minimizing the risk of uterine perforation.*"
- Quote: "*Multi-Modality Imaging Compatibility: Fully compatible with MRI, CT, and X-ray imaging for accurate treatment planning and verification.*"

**Important Facts and Considerations:**

- The device utilizes inflatable ovoids and retractors inflated with sterile saline or water.
- Implant adjustments can be made externally between treatment fractions.

- The system is designed to be anatomically adaptive and conform to patient-specific anatomy.

**Contact Information:**


- Voxelrad@Gmail.com
- Phone: 949-285-3347